

# Guide for digital occlusion

## Interpretation of the OccluSense® data





## INTRODUCTION

Since the beginning of dentistry, the occlusal adjustment is one of the most discussed topics.


Classic articulating papers or foils mark the occlusal contacts in intercuspation. However, the temporal recording of occlusal contacts during dynamic jaw movements or the detailed masticatory pressure distribution are left unrecorded.

OccluSense® allows to visualize and record the entire temporal jaw movement sequence up to the final intercuspation including the relative masticatory pressure distribution of the jaw.

The interpretation of the occlusal pressure, recorded by OccluSense®, is different from classical occlusion test materials as much more information is being provided. This information includes the recording of the occlusal contacts during the slide of the mandible from the initial contact to the maximal intercuspal position.

OccluSense® by Bausch now enables every dentist to record these occlusal situations and evaluate them step by step.

This guide illustrates how OccluSense® recordings can be interpreted to determine an occlusal status for each patient individually.





# CONTENTS




	Page
<b>Occlusion in motion</b> .....	4
<b>OccluSense® technology</b>	Page
Technology: The OccluSense® sensor .....	8
Technology: how the OccluSense® works .....	9
Meaning of colors: Does "red" imply alert? .....	10
Meaning of colors: green contact points .....	11
Meaning of columns: 3-Dimensional visualization of pressure .....	12
Meaning of columns: pressure change ( $\Delta p$ ) visualized by colors .....	13
Additional controls: filtering .....	14
Additional controls: motion controls .....	15
<b>Occlusal situations</b>	Page
Static occlusion .....	16
Dynamic occlusion - Laterotrusion .....	18
Dynamic occlusion - Protrusion .....	20
Splint therapy .....	22




## OCCLUSION IN MOTION

INITIAL CONTACT



SLIDING INTO INTERCUSPATION



0 SEC.

TIME



### MAXIMUM INTERCUSPAL POSITION



20 SEC.

Although occlusion is a process in motion over time, checking the occlusal situation with articulating papers or occlusion test foils is always a snapshot of the final intercuspital position.

Even though classic occlusion test material such as articulating paper or foil mark every occlusal contact, it can not be determined when these contacts have occurred in time. Thus, the occlusal contacts during the temporal mandible movement and the detailed masticatory pressure distribution of the jaw cannot be represented.





## OCCLUSION IN MOTION


INITIAL CONTACT

1



SLIDING INTO INTERCUSPATION

2




0 SEC.

TIME




## OCCUSION IN MOTION

### MAXIMUM INTERCUSPAL POSITION



Occlusense® allows to visualize the entire temporal jaw movement sequence.

Occlusense® records the occlusal contacts during the slide of the mandible from the initial contact to the maximal intercuspal position over the entire course of time.








## TECHNOLOGY


### THE OCCLUSENSE® SENSOR

The OccluSense® sensor is a 60µm thin color coated foil with a printed circuit containing 1018 pressure sensitive pixels which are able to capture 256 levels of pressure. The thin and flexible sensor is able to record occlusal contacts with a low pressure as well as occlusal contacts during dynamic movements of the jaw.

As the sensor is color coated, the exact contact points will be marked on the patient's teeth additionally. OccluSense® is being used like a traditional articulating paper or occlusion test foil but enables the dentist to evaluate the masticatory pressure from the first contact to the maximum intercuspation.



The colors show the change in masticatory pressure of the occlusal contact points in relation to each other and thus visualize the nature of each contact point.



The height of the columns visualizes the relative masticatory pressure between all contact points of the full arch.





## TECHNOLOGY



### HOW THE OCCLUSENSE® WORKS

The Occlusense® recording displays the occlusal situation of the full jaw in different colors. In the 3D view, columns in different heights and colors are shown additionally. The colors are calculated by the relative pressure between the contact points while the height of the columns only shows the masticatory pressure.

The **maximum pressure change ( $\Delta p$ )** between contact points is **red**, the **minimum pressure change ( $\Delta p$ )** between contact points is **green**. The other colors symbolize the values within this range.



$\Delta p$  (Delta  $p$ ) denotes the change ( $\Delta$ ) in pressure ( $p$ ).







# MEANING OF COLORS

## DOES "RED" IMPLY ALERT?

The red color does **not** imply alert. These contact points just show a large pressure change ( $\Delta p$ ) between one contact point and its adjacent contact points.

Standalone red points or columns are typical for:

- highspots
- initial occlusal contacts
- occlusal contacts which are not in intercuspation
- occlusal contacts recorded during dynamic movements






## MEANING OF COLORS

### GREEN CONTACT POINTS

Green contact points show a small pressure change ( $\Delta p$ ) between one contact point and its adjacent contact points.

This is typical for:

- planar contact points
- occlusal contacts on abrasive teeth
- areas of multiple occlusal contact points with similar characteristics






# MEANING OF COLUMNS

## 3-DIMENSIONAL VISUALIZATION OF PRESSURE

The height of the columns visualizes the relative masticatory pressure between all contact points of the full arch.

The visualization of the relative masticatory pressure does not depend on the colors!






# MEANING OF COLUMNS

## PRESSURE CHANGE ( $\Delta p$ ) VISUALIZED BY COLORS

The colors are defined by the pressure change ( $\Delta p$ ) between one contact and all adjacent contacts.






# ADDITIONAL CONTROLS

## FILTERING

The masticatory pressure is being captured in a range between 0 and 256 levels. Every recording stored in the OccluSense®-iPad-App shows a slider for adjusting the pressure level. During the playback of the recordings, this slider can be used by the dentist to eliminate “unwanted” information for a precise evaluation of the contact points. At the same time, **the color information will be recalculated** by taking the remaining contact points into consideration. To avoid the data being altered by electrical noise of the sensor, a threshold is set to a default value of "30 pressure levels".

Technical representation





## ADDITIONAL CONTROLS

### MOTION CONTROLS

Even though occlusion is a movement, articulating papers and foils show a static snapshot. When checking the dynamic occlusion, laterally or protrusive, the markings of articulating papers and occlusion test foils always show the final step of the occlusal movement. With OccluSense®, the motion of the occlusion is being captured from the first contact to the final intercuspation. The occlusal movements are being recorded like a movie, they can be stopped, rewind, forwarded and filtered. Additional viewing options show the percentage of the relative pressure or a quartered pressure force distribution.

The screenshot displays the OccluSense® software interface. On the left is a list of recordings with columns for patient name, frequency, and duration. The selected recording is highlighted in blue. The main area shows a 2D/3D pressure distribution map of a dental arch. The 2D view is a circular map with color-coded pressure zones (green, red, yellow). The 3D view shows a bar chart of pressure distribution. A timeline at the bottom allows for scrubbing through the recording, with markers for current time (02:820) and total time (09:000). Playback controls include 'Rewind one picture', 'Play/Pause', and 'Forward one picture'. A 'Filtering' slider is visible on the right. Red annotations highlight various features: 'Additional viewing options (Percentage, quartered force distribution)', 'Full screen', 'Filtering', 'Current time of recording', 'Scrubbing', 'Total time of recording', and 'Use finger to turn'.

<Back

Diggle, Arnor  
50 Hz 0:05

**Recording**  
**02/14/2020 3:48 PM**

Diggle, Arnor  
100 Hz 0:08

**Recording**  
**02/14/2020 3:45 PM**

Diggle, Arnor  
100 Hz 0:08

**Recording**  
**09/28/2017 3:54 PM**

Diggle, Arnor  
50 Hz 0:09

**Recording**  
**09/28/2017 3:51 PM**

Diggle, Arnor  
50 Hz 0:09

Compare two recordings

**Compare recordings**

2D 3D 2D/3D

Recording 50 Hz 9 seconds

Additional viewing options  
(Percentage, quartered force distribution)

Full screen

Filtering

Use finger to turn

Current time of recording 02:820

Scrubbing

Total time of recording 09:000

Rewind one picture Play/Pause Forward one picture





# OCCLUSAL SITUATIONS

## STATIC OCCLUSION



### Maximum intercuspation

- Regular occlusion on the molar
- Regular occlusion on the premolar
- Premature contact


Occlusion Test Material:  
Bausch Articulating Paper with progressive  
color transfer BK 01 - 200 microns





# OCCLUSAL SITUATIONS

## STATIC OCCLUSION



- Regular occlusion on the molar**  
Contacts with high pressure = high columns higher pressure than adjacent contacts = high  $\Delta p$  -> **RED**
- Regular occlusion on the premolar**  
Contacts with low pressure = low columns no adjacent contacts = high  $\Delta p$  -> **RED**
- Premature contact**  
Contacts with high pressure = high columns, similar  $\Delta p$  between adjacent contacts -> **RED and ORANGE**



# OCCLUSAL SITUATIONS

## DYNAMIC OCCLUSION - LATEROTRUSION

### Laterotrusion right

 Regular occlusion on the molar

 Group guidance

 Premature contact



**Occlusion Test Material:**  
Bausch Progress 100 with progressive color transfer BK 51 - 100 microns blue




**Occlusion Test Material:**  
Bausch OccluSense® Pressure Sensor BK 5025 - 60 microns red





# OCCUSAL SITUATIONS

## DYNAMIC OCCLUSION - LATEROTRUSION



Filter level

- Regular occlusion on the molar**  
Laterotrusion with low pressure = low columns, higher pressure than adjacent contacts = high  $\Delta p$  -> **RED**
- Group guidance**  
Group Guidance with high pressure = high columns, higher pressure than adjacent contacts = high  $\Delta p$  -> **RED**  
Group Guidance with lower pressure = lower columns caused by movement = **RED, ORANGE, GREEN**
- Premature contact**  
Contacts with high pressure = high columns, similar  $\Delta p$  between adjacent contacts -> **RED and ORANGE**




# OCCLUSAL SITUATIONS

## DYNAMIC OCCLUSION - PROTRUSION

### Protrusion

- Anterior canine guidance
- Habitual occlusion -> protrusion
- Premature contact



**Occlusion Test Material:**  
Bausch Arti-Foil Articulating-Foil  
two-sided BK 24 - 8 microns black



**Occlusion Test Material:**  
Bausch OccluSense® Pressure Sensor  
BK 5025 - 60 microns red







# OCCLUSAL SITUATIONS

## DYNAMIC OCCLUSION - PROTRUSION

### Habitual Occlusion



- Anterior canine guidance  
Contacts with low pressure = low columns  
no adjacent contacts = high  $\Delta p$  -> RED
- Habitual occlusion -> protrusion  
Contacts with low pressure = low columns  
no adjacent contacts = high  $\Delta p$  -> RED  
Contact points with lower pressure =  
similar  $\Delta p$  between adjacent contacts ->  
ORANGE, YELLOW, GREEN
- Premature contact  
Contact with high pressure =  
high column, high  $\Delta p$  -> RED




### Protrusion

- Anterior canine guidance  
Contacts with low pressure = low columns  
no adjacent contacts = high  $\Delta p$  -> RED
- Premature contact  
Contact with high pressure =  
high column, high  $\Delta p$  -> RED



# OCCLUSAL SITUATIONS

## SPLINT THERAPY



### Adjustment of splints

○ Static contacts

○ Protrusion movement/Premature contact

○ Laterotrusion movement

○ Protrusion movement

Occlusion Test Material:  
Bausch OccluSense® Pressure Sensor  
BK 5025 - 60 microns red











# OCCUSAL SITUATIONS

## SPLINT THERAPY

Static Occlusion



- **Static contacts**  
Contacts with high pressure = high columns  
(one contact with low pressure = low column)  
no adjacent contacts = high  $\Delta p$  -> **RED**



Dynamic Occlusion (Protrusion)

- **Protrusion movement/Premature contact**  
Contacts with high pressure = high columns  
no adjacent contacts = high  $\Delta p$  -> **RED**
- **Laterotrusion movement**  
Contacts with lower pressure = lower columns  
no adjacent contacts = high  $\Delta p$  -> **RED**
- **Protrusion movement**  
Contacts with lower pressure = lower columns  
no adjacent contacts = high  $\Delta p$  -> **RED**

Supported by:



on the basis of a decision  
by the German Bundestag



[www.occlusense.com](http://www.occlusense.com)

**Bausch**<sup>BK</sup>  
WE MAKE OCCLUSION VISIBLE

Dr. Jean Bausch GmbH & Co. KG | Oskar-Schindler-Straße 4 | D-50769 Köln - Germany | Tel: +49-221-70936-0 | Fax: +49-221-70936-66 | [info@bauschdental.de](mailto:info@bauschdental.de)

Bausch Articulating Papers, Inc. | 12 Murphy Drive, Unit 4 | Nashua, NH 03062, U.S.A. | Tel: +1-603-883-2155 | Tel: 888-6-BAUSCH | Fax: +1-603-883-0606 | [info@bauschdental.com](mailto:info@bauschdental.com)

Bausch Articulating Papers (Australasia) Pty. Ltd | G.P.O. Box 3733, Sydney NSW 2001, Australia | Tel: +61-2-9345-1945 | Fax: +61-2-9345-1955 | [info@bauschdental.com.au](mailto:info@bauschdental.com.au)

Bausch Articulating Papers Japan K. K. | 2nd Floor, 1-4-2, Jonan, Ikedashi | Osaka 563-0025, Japan | Tel: +81 72-737-9501 | Fax: +81 72-737-9502 | [info@bauschdental.jp](mailto:info@bauschdental.jp)

Bausch Importação de Materiais Odontológicos Ltda. | Rua Paulo Eduardo Xavier de Toledo, 379 salas 8 e 9 | 13304-240 Itu-SP, Brasil | Tel: +55 11 3020-9263 | [vendas@bauschbrasil.com.br](mailto: vendas@bauschbrasil.com.br)

Bausch Articulating Papers Co., Ltd. | #1004, The O Ville B/D | 344, Jong-Ro, Jongno-Gu | Seoul 03114 | Korea | [info@bauschdental.kr](mailto:info@bauschdental.kr) | [www.bauschdental.kr](http://www.bauschdental.kr)

Bausch and OccluSense® are trademarks of Dr. Jean Bausch GmbH & Co. KG, registered in the EU and other countries. Apple, iPad, iPad Air are trademarks of Apple Inc., registered in the U.S. and other countries.