Prior to use, carefully read the instructions for use. EN

# GC Fuji TRIAGE® CAPSULE

RADIOPAQUE GLASS IONOMER PROTECTION AND TEMPORARY RESTORATIVE MATERIAL IN CAPSULES

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

1. Fissure protection.
2. Root surface protection.
3. Hypersensitivity prevention and control.
4. Protection for immature enamel.
5. Temporary filling of endodontic access.
6. Intermediate restorations.

### CONTRAINDICATIONS

1. Pulp capping
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

PINK Shade is command set, WHITE Shade is chemical set only .

Powder / Liquid Ratio (g / g)	0.30 / 0.15
Mixing Time (sec.)	10"
Working Time (min., sec.)	140"
Net Setting Time (min., sec.)	230"
Final Finishing Commencing Time	600"
Final Finishing Commencing Time if light cured (PINK Command set)	400"

Test conditions : Temperature (23±1°C), Relative humidity (50±5%)  
 ISO 9917-1 : 2003 (Dental water-based cements) (Restorative cements)

### A. FISSURE PROTECTION ROOT SURFACE PROTECTION HYPERSENSITIVITY PREVENTION AND CONTROL PROTECTION FOR IMMATURE ENAMEL

1. Preparation of the tooth surfaces (e.g. for fissure protection or root surface protection)

a) After cleaning the tooth surfaces (prophylaxis with pumice and water) in usual manner, rinse thoroughly with water. Avoid aggravating the operculum.  
 Note :  
 If extra retention is desired, application of GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) is recommended.

b) Dry by blotting with a cotton pellet or gently blowing with air syringe (Fig. A-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).

### 2. Mixing

a) Before activation, shake the capsule or tap its side 2 or 3 times on a hard surface to loosen the powder (Fig. A-2).  
 b) To activate the capsule, push the plunger until it is flush with the main body (Fig. A-3).  
 c) Immediately place the capsule into a GC Capsule Applier and click the lever once (Fig. A-4). The capsule is now activated.

Note :  
 The capsule should be activated just before mixing and used immediately.  
 d) Immediately remove the capsule from the applier and set the capsule into a capsule mixer or amalgamator. Mix for 10 seconds at high speed (approximately 4,000 RPM) (Fig. A-5).

3. Placement  
 a) Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier.  
 b) Make two clicks to prime the capsule then syringe (Fig. A-6). The working time is 1 minute 40 seconds from the start of mixing at 23°(73.4°F). Higher temperatures will shorten working time.

c) Extrude the mixture onto the tooth surface and use a brush to spread a thin film of GC Fuji TRIAGE® directly over the root surface or hypersensitive area or over the occlusal surface and into the pits and fissures.  
 Note :  
 1) If a faster set is desired, use a visible light curing device\* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.

2) To adjust the direction of the nozzle, hold the applier with the capsule towards you and turn the capsule body.  
 3) To remove the use capsule, push the applier release button. Twist the capsule and pull upwards.

d) After placement, when the material starts to lose the glossy appearance (or after curing with the light curing device), apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) to the sealed area and the margins using a cotton pellet or sponge (Fig. A-8).  
 e) Finishing under air water spray can be performed 6 minutes from start of mix (chemically set) or 4 minutes if light cured. Use a superfine diamond bur or a silicone finishing point (Fig. A-9).

f) Apply GC Fuji VARNISH or GC Fuji COAT LC to the area again (Fig. A-10).

### B. TEMPORARY FILLING OF ENDODONTIC ACCESS

1. Cleaning the pulp chamber  
 a) After appropriate pulp treatment, clean and gently dry the pulp chamber with an air syringe (Fig. B-1).  
 b) Fill the chamber with a cotton pellet (Fig. B-2).

### 2. Mixing

See directions in A (above), Section 2.  
 3. Placement  
 a) See directions in A, Section 3, a and b.  
 b) Extrude the mixture directly over the cotton pellet (Fig. B-3).  
 Note :  
 If a faster set is desired, use a visible light curing device\* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.

c) Moisture protection and Finishing  
 See directions in A, Section 3, d, e and f.

### C. INTERMEDIATE RESTORATIONS

1. Caries removal  
 a) Remove caries and loose debris with hand instruments.  
 b) For better retention, it is recommended to gently clean the carious surface with GC CAVITY CONDITIONER for 10 seconds or GC DENTIN CONDITIONER for 20 seconds.

c) Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. C-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).

### 2. Mixing

See directions in A (above), Section 2.  
 3. Placement  
 a) See directions in A, Section 3, a and b.  
 b) Extrude the mixture directly over the prepared lesion or tooth surface (Fig. C-2). Form the contour and if possible cover with a matrix.  
 Note :  
 If a faster set is desired, use a visible light curing device\* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.

c) See directions in A, Section 3, d, e and f.

\* NOTE :  
 The initial set of GC Fuji TRIAGE® PINK can be accelerated using the energy from a dental halogen light curing device.

### SHADE

Pink, White  
**STORAGE**  
 Store at temperature of 4-25° (39.2-77.0°F)  
 (Shelf life : 2 years from date of manufacture)

**PACKAGE**  
 50 capsules  
 Average contents per capsule : 0.30g powder and 0.15g (0.12mL) liquid  
 Minimum net volume of mixed cement per capsule : 0.13mL

**CAUTION**  
 1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.  
 2. In case of contact with eyes, flush immediately with water and seek medical attention.

US Patent : 6264472 6756421  
 UK Patent : 2353042 2357773  
 France Patent : 2797396 2799954  
 Australia Patent : 768901 775349

Last revised : 07/2008

MANUFACTURED BY  
 GC CORPORATION  
 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY  
 GC CORPORATION  
 76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.  
 Research-Park, Interleuvenlaan 13, B-3001 Leuven, Belgium  
 TEL: +32. 16. 39. 80. 50.

GC AMERICA INC.  
 3737 West 127th Street, Alsip, IL 60803 U.S.A.  
 TEL: +1. 708. 977. 0800

GC ASIA DENTAL PTE. LTD.  
 19 Loyang Way, #06-27 Singapore 508724  
 TEL: +65 6546 7588

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EN Prior to use, carefully read the instructions for use.

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### RECOMMENDED INDICATIONS

1. Fissure protection.
2. Root surface protection.
3. Hypersensitivity prevention and control.
4. Protection for temporary restorations.
5. Temporary filling of endodontic access.
6. Intermediate restorations.

### CONTRAINDICATIONS

1. Patients with severe allergic reactions to the product and/or any of its components.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and consult a physician.

### DIRECTIONS FOR USE

PINK Shade is command set, WHITE Shade is chemical set only.

Powder / Liquid Ratio (g / g)	0.30 / 0.15
Mixing Time (min., sec.)	10*
Working Time (min., sec.)	140*
Net Setting Time (min., sec.)	230*
Final Finishing Compaction Time	6'00"
Final Finishing Compaction Time if light cured (PINK Command set)	4'00"

Test conditions: Temperature (23±1 °C), Relative humidity (50±5%)

### ISO 9917-1: 2003 (Dental water-based cements) (Restorative cements)

### A. FISSURE PROTECTION

#### ROOT SURFACE PROTECTION HYPERSENSITIVITY PREVENTION AND CONTROL PROTECTANT FOR IMMATURE ENAMEL

1. Preparation of the tooth surfaces (e.g. for fissure protection or root surface protection)
- a) After cleaning the tooth surfaces (prophylaxis with fluoride and water) in usual manner, rinse thoroughly with water. Avoid aggravating the operculum.

If extra retention is desired, application of GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) is recommended.

- b) Dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. A-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).

2. Mixing
- a) Before activation, shake the capsule or put its side or 3 times on a hard surface to loosen the powder.

- b) To activate the capsule, push the plunger until it is flush with the main body (Fig. A-3).

- c) Immediately place the capsule into a GC Capsule Applicator and click the lever once (Fig. A-4). The capsule is now activated.

The capsule should be activated just before mixing and used immediately.

- d) Immediately remove the capsule from the applicator and insert it into a capsule mixer or amalgamator.

- e) Run for 10 seconds at high speed (approximately 4,000 RPM) (Fig. A-5).

The capsule should be activated just before mixing and used immediately.

- f) Immediately remove the capsule from the applicator and insert it into a capsule mixer or amalgamator.

- g) Run for 10 seconds at high speed (approximately 4,000 RPM) (Fig. A-5).

3. Placement
- a) Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applicator.

- b) Make two clicks to prime the capsule then syringe (Fig. A-6). The working time is 1 minute 40 seconds from the start of mixing at 23 (73.4 °F). Higher temperatures will shorten working time.

- c) Push the mixture directly onto the tooth (Fig. A-7) then use a brush to spread a thin film of GC Fuji TRIAGE® directly over the root surface or over the occlusal and incisal surfaces and into the pits and fissures.

Note: If a master set is desired, use a visible light curing device for 20-40 seconds. Place light source as close as possible to the cement surface. This method is recommended only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.

- d) To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- e) To remove the used capsule, push the applicator release button. Twist the capsule, pull up upwards.

- f) After placement, when the material starts to lose its shine, a thin film of curing device (or a light curing device), apply GC Fuji VARNISH (blow dry the GC Fuji COAT LC light cure) to the sealed areas and the margins using an water pellet or sponge (Fig. A-8).

- g) Finishing under air water spray can be performed 6-8 minutes after placement on the preparation.

- h) Use a superfine diamond bur or a silicone finishing point. Use GC Fuji COAT LC to the area again (Fig. A-10).

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1. Cleaning the pulp chamber

- a) After appropriate pulp treatment, clean and gently dry the pulp chamber with an air syringe (Fig. B-1).

- b) Fill the chamber with a cotton pellet (Fig. B-2).

### 2. Mixing

- a) See directions in A (above), Section 2.

### 3. Placement

- a) See directions in A, Section 3, a and b.

- b) Extrude the mixture directly over the cotton pellet (Fig. B-3).

### 4. Finishing

- a) See directions in A, Section 3, d, e, and f.

### C. INTERMEDIATE RESTORATIONS

1. Caries removal

- a) Remove caries and loose debris with hand instruments.

- b) For better retention, it is recommended to gently clean the carious surface with GC CAVITY CONDITIONER for 10 seconds or GC DENTIN CONDITIONER for 20 seconds.

- c) Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. C-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).

### 2. Mixing

- a) See directions in A (above), Section 2.

### 3. Placement

- a) See directions in A, Section 3, a and b.

- b) Extrude the mixture directly over the prepared lesion or tooth surface (Fig. C-2). Form the contour and if possible cover with a matrix.

Note: If a faster set is desired, use a visible light curing device\* for 20-40 seconds. Place light source as close as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.

- c) Moisture protection and finishing

- a) See directions in A, Section 3, d, e, and f.

### \* NOTE:

The initial set of GC Fuji TRIAGE® PINK can be accelerated using the energy of a dental halogen light curing device.

### SHADE

Pink, White

### STORAGE

Store at a temperature of 4-25° (39-77.0°F)

(Shelf life: 2 years from date of manufacture)

### PACKAGE

50 capsules

Contents per capsule: 0.30g powder and 0.15g (0.12mL) liquid

Minimum net volume of mixed cement per capsule: 0.13mL

### CAUTION

1. In case of contact with oral tissue or skin, remove immediately with water or cotton soaked in alcohol. Flush with water.

2. In case of contact with eyes, flush immediately with water and seek medical attention.

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