

Glass Ionomers Restoratives and Liners

Vitremer^a

An advanced glass ionomer system that offers a three-way cure to meet your needs for core buildup and restorative applications.

Indications

- Class III and Class V restorations.
- Root Caries.
- Pediatric I or II.
- Base/liner.
- Class II laminate or open sandwich.
- Core buildup.

Product Benefits

- Versatile system can be light-cured or self-cured to meet your application needs.
- Unique dark-cure chemistry allows for bulk placement, which eliminates the need for layering and saves time.
- Displays a true glass ionomer cure for long-term fluoride release and adhesion to tooth structure, which provides increased doctor confidence.
- As a restorative material, it provides good handling properties and extended working time in the mouth for easier use.
- Available in a variety of shades to better match existing tooth structure for esthetic restorations. Pedito shade is available for pedodontic applications and Blue shade for core buildup applications.



Product Information	Item Number
Introductory Kit:3303I
5g Powder each of A3, A4, C2, C4, PEDO and BLUE; 2-D8ml Liquid; 6.5ml Primer; 6.5ml Finishing Gloss; Shade Guide and Accessories	
Trial Kit:3303MPA3
5g Powder of A3 Shade; 2.5ml Liquid; 2ml Primer; 2ml Finishing Gloss, Delivery Tips with Pistons (10 each); Spoon and Mixing Pad	
Trial Kit:3303MPB
5g Powder of Blue Shade; 2.5ml Liquid; 2ml Primer; Delivery Tips with Pistons (10 each); Spoon and Mixing Pad	
Trial Kit:3303MPP
5g Powder of Pedito Shade; 2.5ml Liquid; 2ml Primer; 2ml Finishing Gloss; Delivery Tips with Pistons (10 each); Spoon and Mixing Pad	

Product Information	Item Number
Refills:	
1-D 9gA3 Shade Powder, Spoon3303A3
1-D 9gA3.5 Shade Powder, Spoon3303A3.5
1-D 9gA4 Shade Powder, Spoon3303A4
1-D 9g B2 Shade Powder, Spoon3303B2
1-D 9g B3 Shade Powder, Spoon3303B3
1-D 9g C2 Shade Powder, Spoon3303C2
1-D 9g C4 Shade Powder, Spoon3303C4
1-D 9g PEDO Shade Powder, Spoon3303PEDO
1-D 9g BLUE Shade Powder, Spoon3303BLUE
1-D 8ml Liquid3303L
1-D 6.5ml Primer3303P
1-D 6.5ml Finishing Gloss3303FG
Delivery Tips with Pistons, 50 each3303DT